

# Specifications for Agarose Gels

CATALOG NO.	DESCRIPTION	GEL SIZE	BUFFER VOLUME	RECOMMENDED VOLTAGE	APPROXIMATE RUN TIME
3486-5000	DNA Plus <sup>®</sup> mini gel	7 x 8 cm	400 ml	100 V	30-60 minutes
3487-0000	DNA Plus <sup>®</sup> mini gel	9 x 11 cm	600 ml	125 V	60-90 minutes
3487-5000	DNA Plus <sup>®</sup> mini gel	12 x 14 cm	800 ml	150 V	60-120 minutes
3487-3000	DNA Plus <sup>®</sup> , self-circulating	12 x 14 cm	1000 ml	150 V	60-240 minutes
3431-4000	DNA Plus <sup>®</sup> , high capacity	23 x 14 cm	800 ml	150 V	30-60 minutes
3482-0000	DNA Plus <sup>®</sup> , large	20 x 25 cm	2300 ml	150-250 V	1-4 hours

