

ErgoOne® & TipOne® Compatibility Chart

ErgoOne® Volume Range	TipOne® Pipet Tips	Sterile Filter Tips	Filter Tip Refills	Racks	Refill Wafers	Stacks	Starter Systems	Sterile Racks	Bulk Tips
0.1-2.5 µl 0.5-10 µl ¹ 2-20 µl micro tip cone ²	10 µl	1121-3810	1121-2710	1111-3800	1111-3700	1111-3200 1111-3210	1111-3500	1111-3810	1111-3000
		1181-3810	1181-3710	1161-3800	1161-3700	--	1161-3500	--	--
	10 µl XL	1120-3810	1120-3710	1110-3800	1110-3700	--	1110-3500	1110-3810	1110-3000
		1180-3810	1180-3710	1160-3800	1160-3700	--	1160-3500	--	--
2-20 µl standard tip cone 10-100 µl 20-200 µl	20 µl beveled	1120-1810	1120-1710	--	--	--	--	--	--
		1180-1810	1180-1710	--	--	--	--	--	--
	20 µl Profile Point™	1123-1810	1123-1710	--	--	--	--	--	--
		1183-1810	1183-1710	--	--	--	--	--	--
	100 µl beveled	1120-1840	1120-1740	--	--	--	--	--	--
		1180-1840	1180-1740	--	--	--	--	--	--
	100 µl Profile Point™	1123-1840	1123-1740	--	--	--	--	--	--
		1183-1840	1183-1740	--	--	--	--	--	--
	200 µl graduated	1120-8810	1120-8710	--	--	--	--	--	--
	200 µl beveled	--	--	1111-0800	1111-0700	1111-0200	1111-0500	1111-0810	1111-0000
	200 µl yellow	--	--	1111-0806	1111-0706	1111-0206	1111-0506	1111-0816	1111-0006
	200 µl beveled	--	--	1111-1800	1111-1700	1111-1200 1111-1210	1111-1500	1111-1810	1111-1000
		1180-8810	1180-8710	1161-1800	1161-1700	--	1161-1520	--	--
	200 µl yellow, beveled	--	--	1111-1806	1111-1706	1111-1206	1111-1506	1111-1816	1111-1006
200 µl Profile Point™	--	--	1110-1800	1110-1700	1110-1200	1110-1500	1110-1810	1110-1000	
30-300 µl ³	300 µl graduated	1120-9810	1120-9710	1110-9800	1110-9700	--	1110-9500	1110-9810	1110-9000
		1180-9810	1180-9710	1160-9800	1160-9700	--	1160-9500	--	--
100-1000 µl	1000 µl	1126-7810	1126-7710	--	--	--	--	--	--
	1000 µl graduated	--	--	1111-2820	1111-2720	--	1111-2520	1111-2830	1111-2020
	1000 µl blue, graduated	--	--	1111-2821	1111-2721	--	1111-2521	1111-2831	1111-2021
	1250 µl XL	1122-1830	1122-1730	1112-1820	1112-1720	--	1112-1500	1112-1830	1112-1020
1182-1830		1182-1730	1161-1820	1161-1720	--	1161-1520	--	--	
0.5-5 ml ⁴	5 ml graduated	1050-0810	--	1050-0700	--	--	--	1050-0710	1050-0000

Original TipOne pipet tips

TipOne RPT ultra low retention tips

¹Tips in these rows are recommended for 0.5-10 µl single channels

²2-20 µl micro tip cone pipet is designed to use ultra micro tips only. It will accept any of the ultra micro tips above.

³30-300 µl will also accept 200 µl tips.

⁴See page A23 for more information on 5 ml tips.