




# x-tracta™ Disposable Gel Extraction Tool

- Eliminates cross contamination
- Minimizes UV exposure and damage
- Cuts uniform gel slices
- Replaces hazardous blades

The x-tracta agarose gel tool gives you cleaner DNA in less time and with less effort. This disposable device cleanly cuts and holds uniform gel slices, eliminating cross-contamination from razor blades, tweezers, or other tools. The process is quick, reducing damaging UV exposure. Gel bands can be released directly into microcentrifuge tubes for the next extraction step.



## INSTRUCTIONS FOR USE



In accordance with good laboratory practice, use the x-tracta tool once, then discard it to prevent cross-contamination.  
 Note: Do not squeeze the bulb until step 4. At step 4, use a quick snap-action squeeze to release the gel.