

96-WELL PCR PLATE SELECTION GUIDE

**Standard
Well
Depth**


Skirt	Semi	Semi	Semi	None	None
Additional Features	Contoured sides, rimmed wells	Raised contoured sides, rimmed wells	Straight sides, rimmed wells	Rimmed wells	Chimney wells
Notched Corner	A12	A12	A12	A12	H12
Working Volume (per well)	1-240 μ l	1-240 μ l	1-240 μ l	1-240 μ l	1-240 μ l
Maximum Volume	300 μ l	300 μ l	300 μ l	300 μ l	300 μ l
Color Options	Natural, Assorted, White, Black	Natural, Assorted	Natural	Natural, Assorted, White	Natural, Assorted
High Contrast Alphanumerics	Yes	Yes	Yes	Yes	Yes
Natural Plate #	1402-9100	1402-9300	1402-9200	1402-9596	1402-8600
Catalog Page #	CI	CI	C6	C9	C8

Are you working with very low volumes? Try our low profile PCR tubes or plates.

**Low
Profile
Wells**


Skirt	Semi	Semi	Semi	None	Full
Additional Features	Contoured sides, rimmed wells	Raised straight sides, rimmed wells	Straight sides, rimmed wells	Rimmed wells	Rimmed wells
Notched Corner	A12	A1	H12	H12	Natural A12 White H1
Working Volume (per well)	1-140 μ l	1-140 μ l	1-140 μ l	1-140 μ l	1-140 μ l
Maximum Volume	200 μ l	200 μ l	200 μ l	200 μ l	200 μ l
Color Options	Natural	Natural, White	Natural, White	Natural, White	Natural, White
High Contrast Alphanumerics	Yes	Yes	Yes	Yes	Yes
Natural Plate #	1402-8400	1402-9100	1402-8900	1402-9500	1402-9800
Catalog Page #	CI	C6	C6	C9	C10